

Follow these steps to install the new version.

1. Perform a backup. It can either be an internal backup or a server backup. You don't have to perform a CD backup so make sure the "Copy To Media" option is un-checked.
 2. Make sure everyone has exited the program.
 3. Go to the Revisions screen on our website (www.drmnet.com). Shift and click the link to download the updates (remember to hold your shift key down while clicking).
-
- a. Login name is: drmwuser
 - b. Password: (You will not see characters as you type.) Use this month's password.
 - c. After you click the Login button you will see a list of files.
 - d. If you see a folder named drmwuser then double-click on it otherwise you are already in the right place.
4. Perform the Update.
 - a. Double-click the file `drm_180_596.exe` and choose to run it.
 - b. After a moment the DRM Workstation Update Wizard will run. If you see the following message, "This will update all these companies to the current version", Click <Proceed>
 5. You need to perform the update on EVERY computer where DRM is used.